

## Control cell lysates Hmgb1 -/-

**Product Number:** (depends on batch)  
**Expiration date:** (depends on batch)  
**Batch number:** (depends on batch)  
**Batch concentration:** (depends on batch)


### **Product Description:**

Cell lysates from wild type or HMGB1-deficient cells are essential to verify the specificity of antibodies in Western blots, or to test the efficacy of siRNAs directed against HMGB1.

We provide aliquots containing 50 microgram each of total cell protein, from cells lysed directly in SDS-PAGE loading buffer. The aliquots are lyophilized, and shipped at room temperature.

**Storage:** 2-8°C.

**This product is for research use only**



**Fig. 1.** Western Blot